

Data Sheet THERMARMOUR



Emergency Blanket


Product Code	OB000150200---TAEM55
GTIN/UDI	5061035931467
Colour(s)	Emergency Blue Outer with White Inner
Size	150cm x 200cm (can be cut to size if required)
Material	Non-conductive bilaminate foil & needle-punched polyester fleece
Care	Patient Specific, cannot be laundered
Buy Online	www.interweavetextiles.com/products/patient-warming/

Designed with the goal of maintaining and restoring normothermic temperatures in emergencies, THERMARMOUR is suitable for patients at risk of hypothermia. Other scenarios include, but are not limited to pre-warming for surgeries and temperature management of trauma casualties.

Key Features	<ul style="list-style-type: none"> • Latex Free, Recyclable & X-ray Compatible. • Easily transportable standalone device. • No waiting to use once opened; a bonus for first responders. • A cost and time efficient passive warming method.
---------------------	--

Usage

Open the vacuumed package and THERMARMOUR Blanket.
Unfold the blanket fully, and if required, cut to the desired size.
Wrap around the distressed party with the foil side outwards.



Packing

Manufactured in our UK facility, THERMAMOUR is a readily available, sustainable solution.

Per Carton	30, Individually Vacuum Packed
Carton Size	58 x 40 x 45.5cm
Pallet Quantity	600

Address: Interweave House, Old Power Way, Lowfields Business Park, Elland, HX5 9DE
Tel: 01422 411463 | sales@thermarmour.com | www.thermarmour.com

THERMARMOUR is a trading brand of Interweave Textiles Ltd. Interweave Textiles Ltd is a company registered in England & Wales: 05267141


interweave
healthcare



Data Sheet Storing and Handling

Keep dry and store out of direct sunlight until needed. THERMARMOUR products are for single patient use and should be discarded once soiled. If there is any damage identified upon receipt, please contact sales@thermarmour.com before use.

Technical Drawing



Risks and Precautions

Following the trial and testing period, no known residual risks or side effects have been identified. Please ensure all components including packaging are recycled or disposed of safely, as misuse may cause suffocation.

If there are any safety concerns with this device, please inform THERMARMOUR and MHRA.

Address: Interweave House, Old Power Way, Lowfields Business Park, Elland, HX5 9DE
Tel: 01422 411463 | sales@thermarmour.com | www.thermarmour.com

THERMARMOUR is a trading brand of Interweave Textiles Ltd. Interweave Textiles Ltd is a company registered in England & Wales: 05267141

