

Bed Ladder Strap

Material: Polyester

Colour: Blue

Size: 1.9m long

Product codes: PH801

Washable: Not suitable for laundering. Patient specific use.

Usage: The ladder strap ties onto the bed leg or a bed rail at the

foot of the bed and can be pinned with safety pins to the top layer of the bedding to keep it in reaching distance overnight. There are 5 handholds that act as ladder rungs.



Interweave Textiles Ltd, Interweave House, Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333