

Patient Specific Repositioning Gloves

Material: Antibacterial, antifungal 65 gsm nylon with silicone coating

on the outside and anti-slip inside.

Colour: Red

Size: 75cm x 20cm

Product code: PH803-LGE---REDN07

Washable: Not washable. Disposable.

Usage: Used for a variety of small turning and positioning moves.

Can be used by a carer to move a patient across a bed, can also be used to smooth bed sheets under a patient. Makes movement smoother and easier for both patient and carer.



Interweave Textiles Ltd, Interweave House, Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333