

Handling Slings

- Material:** Polyester outer cover, polyester foam inner and polyester webbing straps.
- Product codes:** SL005-22-62---BLUV24 Box of 10
- Washable:** Launderable at thermal disinfection temperatures, if required. Do not tumble dry. The lining may be wiped with a suitable disinfectant.
- Usage:** Carers use this sling to assist a patient in moving from a sitting to a standing position, by placing the sling around the patient's lower back. There are two handholds to give the carer better grip. The sling should only be used in this way if the patient can weight-bear. Otherwise this sling is useful for performing leg transfers by firstly placing the sling under the patient's lower leg, then moving.

