

Airpal Air Assisted Transfer Equipment

Patient specific lateral transfer mattress

Material: Advanced fabrics with a coating that can be wiped clean with standard antibacterial cleaner.

Size: 198cm long by 86cm wide

Product code: AP020-86198

Washable: Can be laundered to last throughout a patient's stay in hospital.

Usage: The mattress is inflated by the Airpal air supply unit, providing a comfortable lateral transfer for patients with this air assisted lateral patient transfer system. Compatible with XRay and MRI. Snaps are nickel plated brass. Safe working load 1,000kg (or 157 stones).



Re-usable patient lateral transfer mattress

- Material:** Advanced fabrics with a coating that allows for laundering.
- Size:** 198cm long x 86cm wide
- Product code:** AP000-86198
- Washable:** Can be laundered and re-used for many patients.
- Usage:** The mattress is inflated by the Airpal air supply unit, providing a comfortable lateral transfer for patients with this air assisted lateral patient transfer system. Compatible with XRay and MRI. Snaps are nickel plated brass. Safe working load 1,000kg (or 157 stones).



Airpal air supply stand

Material: Manufactured from hospital-grade components, the main stand pole is stainless steel. The base is constructed from corrosion resistant polished cast metal. Provided with two air supply holders and a dishwasher-friendly plastic storage basket.

Height: 116cm

Footprint: 61cm x 61cm (space taken up by the casters on the floor)

Basket size: 31cm wide x 19cm across x 25cm deep

Product code: AP300

Usage: The Airpal stand is for use with the Airpal patient transfer system. Store the air supply in the stand and use the storage basket for the other accessories. Also includes an IV pole.



Airpal air supply

Material: The filter canister is made from non-corrosive stainless steel alloy, and the motor enclosure is high impact moulded plastic. Provided with integrated cable holder and a bed/stretcher hook. Supplied with 180cm long flexible PVC hose. User manual supplied with the air supply.

Height: 28cm

Footprint: 23cm x 23cm (space taken up by the air supply on the floor)

Product code: AP300

Usage: The Airpal air supply hose is connected to the mattress, clip it in place with the poppers then switch on the air flow by pressing the large button on the air supply machine.



Interweave Textiles Ltd, Interweave House, Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333

Airpal mattress cover (sani-liner)

Material: Advanced fabrics with a coating that allows for laundering.

Size: 198cm long x 86cm wide

Product code: AP100

Usage: The Airpal sani-liner attaches to the mattress with non-ferrous poppers. All materials are MRI-friendly. The sani-liner fully covers the mattress when inflated. The sani-liners provide a quick turn-around between patients, reduced laundry costs and an extended life for the re-usable mattress.



Interweave Textiles Ltd, Interweave House, Old Power Way, Lowfields Business Park, Elland, West Yorkshire, UK, HX5 9DE Tel: 01422 372333

Document valid as at 29/07/2019

Page 5 of 5 - Revision no. 5