

Stocking Slider

- Material:** PU coated polyester
- Colour:** Red
- Size:** One size – 22cm x 70cm
- Product code:** PH802-22-70---REDA06 (NHS Supply Chain code: VTS117)
- Washable:** Not launderable – patient specific
- Usage:** Used to fit tight-fitting, open toe compression stockings on patients. Slide the stocking slider onto the foot with the handle end at the toes, then slide the compression stocking over the top and onto the leg. Once the compression stocking is in place, remove the stocking slider by pulling the black handle. Also assists with fitting compression sleeves onto an arm.

